

# Micro-Flask by Duetz, cultivation in microtiter plates

The Micro-Flask system facilitates reproducible and reliable culturing on microtiter plates. The system consists of sandwich covers, cover clamps and cryo-replicator. The Micro-Flask enables a single person to grow and test thousands of strains simultaneously with a minimum of repetitive handling.




## Features













- Conversion of 24 and 96 microtiter plates (both deep- and low-well plates) into individual micro-reactors
- Low and uniform evaporation rates for every well
- Sterile barrier for individual wells prevents cross contamination
- Oxygen transfer rates similar to shake flasks in standard orbital shakers
- Simultaneous and reproducible sampling of 96 frozen glycerol stocks

## Applications

- High throughput screening and distribution of mutant and construct libraries e.g. in *E. coli* or yeast
- Metabolic flux studies and high-throughput screening for high activity prokaryotic or eukaryotic mutants
- Comparative studies, e.g. clinical isolates
- Growth medium optimization for cell lines or production strains



# Specifications

Type of microtiter plate	Well volume	Culture volume	Orbital shaking frequency	Shaking amplitude	O <sub>2</sub> -transfer rate (30°C, air, 1 bar)	Headspace refreshment rate	Evaporation rate per well (at 30°C)	Mixing pattern at 300 rpm	
								ampl.	ampl.
24-square deep-well polypropylene, 17x17 mm, depth 40 mm		2500 µl	300 rpm	50 mm	51 mmol O <sub>2</sub> / l / h	2.5 ml / min	<b>50% humidity:</b> 50 µl H <sub>2</sub> O per day <b>75% humidity:</b> 25 µl H <sub>2</sub> O per day		
		2500 µl	300 rpm	25 mm	39 mmol O <sub>2</sub> / l / h	(1 VVM)		2500 µl	2500 µl
		2500 µl	220 rpm	50 mm	35 mmol O <sub>2</sub> / l / h				
		4000 µl	300 rpm	50 mm	24 mmol O <sub>2</sub> / l / h	2.5 ml / min			
		4000 µl	220 rpm	25 mm	24 mmol O <sub>2</sub> / l / h	(0.6 VVM)			
24-round low-well polystyrene, Ø 16 mm, depth 18 mm		750 µl	300 rpm	50 mm	40 mmol O <sub>2</sub> / l / h	1.1 ml / min	<b>50% humidity:</b> 30 µl H <sub>2</sub> O per day <b>75% humidity:</b> 15 µl H <sub>2</sub> O per day		
		750 µl	300 rpm	25 mm	25 mmol O <sub>2</sub> / l / h	(1.4 VVM)		1000 µl	1000 µl
		1000 µl	300 rpm	50 mm	30 mmol O <sub>2</sub> / l / h	1.1 ml / min			
		1000 µl	300 rpm	25 mm	19 mmol O <sub>2</sub> / l / h	(1.1 VVM)			
96-square deep-well polypropylene, 8x8 mm, depth 40 mm		500 µl	300 rpm	50 mm	38 mmol O <sub>2</sub> / l / h	1 ml / min	<b>50% humidity:</b> 22 µl H <sub>2</sub> O per day <b>75% humidity:</b> 11 µl H <sub>2</sub> O per day		
		500 µl	300 rpm	25 mm	12 mmol O <sub>2</sub> / l / h	(2 VVM)		750 µl	750 µl
		750 µl	300 rpm	50 mm	24 mmol O <sub>2</sub> / l / h	1 ml / min			
		750 µl	300 rpm	25 mm	7 mmol O <sub>2</sub> / l / h	(1.3 VVM)			
		1000 µl	300 rpm	50 mm	18 mmol O <sub>2</sub> / l / h	1 ml / min			
1000 µl	300 rpm	25 mm	3 mmol O <sub>2</sub> / l / h	(1 VVM)					
96-round low-well polystyrene, Ø 6.5 mm, depth 11 mm		100 µl	300 rpm	50 mm	39 mmol O <sub>2</sub> / l / h	250 µl / min	<b>50% humidity:</b> 6 µl H <sub>2</sub> O per day <b>75% humidity:</b> 3 µl H <sub>2</sub> O per day		
		100 µl	300 rpm	25 mm	20 mmol O <sub>2</sub> / l / h	(2.5 VVM)		150 µl	150 µl
		150 µl	300 rpm	50 mm	32 mmol O <sub>2</sub> / l / h	250 µl / min			
		150 µl	300 rpm	25 mm	16 mmol O <sub>2</sub> / l / h	(1.7 VVM)			
		200 µl	220 rpm	50 mm	12 mmol O <sub>2</sub> / l / h	250 µl / min			
200 µl	300 rpm	25 mm	12 mmol O <sub>2</sub> / l / h	(1.3 VVM)					

