

PalPharMassage™

The PalPharMassage has been developed to assist with the process infeed or tipping of powdered or granular product in drums which has hardened or settled in storage and transit.



Mobile PalPharMassage drum conditioning system.



PalPharMassage cylinders can be adjustable to suit different drum heights and powder fills within the drum.



Cylinders can be easily adjusted using stainless steel hand wheels. There is also an option to adjust the backstop for drums of different diameters.

Product Description

The PalPharMassage has been designed to be fully flexible around drums of all sizes, plastic fibre and metal. The drum sits on a rotary table so after each cycle the drum can be rotated. An adjustable back stop allows penetration depth of the cylinder into the drum to be adjusted to suit

the material requirement.

The PalPharMassage can be a simple mobile product conditioning solution or integrated into a turn key fully contained drum discharge solution with PLC control.

Features and Benefits

- Requires only a clean dry compressed air supply (unless PLC control is required)
- Simple to use
- Adjustable cylinder height and backstop options to suit a range of drums
- Mobile or Fixed Options
- cGMP construction with 304 or 316 options available
- Can be used with PharmaVac/
- PharmaTip to offer turn key manual handling and drum tipping solution
- Heavy Duty Construction
- Safety Guarding Options
- CE Label

Product:

Drums, Kegs, Bottles, Containers

Capacity

Up to 300kg

Applications

Conditioning product to assist with the process infeed or tipping of powdered or granular product in drums which has hardened or settled in storage and transit

Finishes

GMP construction with 304 or 316 options available

Atex options

Atex options available for zone 1/21 2/22 T3,T4,T5 Class 1 Div 1/Group C + D

Documentation

FAT and IQ/OQ documentation