


***Romaco to become number one in the Powder Filling systems***

Totally focused on becoming the number one in Powder Filling Systems, Romaco has launched a new, innovative range of powder filling machine.

In order to meet new market requirements, Romaco is expanding its machinery portfolio to include the new MicroMaxX.

The new concept behind the technical design of the MicroMaxX Series leads Romaco closest than ever to the number one position in the world of aseptic powder process.



Technical Data	MicroMaxX 18		MicroMaxX 24	
Vials Diam., mm (min./max)	16 - 52	16 - 32	16 - 52	16 - 32
Vials H., mm (min./max)	35 - 110	35 - 75	35 - 110	35 - 75
Max. output (vph) Single dose per dosing disc, no C.W.	9.000	18.000	12.000	24.000
Max output (vph) Single dose per dosing disc, 100% C.W.	6.000	12.000	12.000	24.000
Min./Max dose (single shot per dosing disc), mg Min./Max dosage can change in relation to powder specific weight	30 – 1.500			
Dosing disc	1 / 2			
Dosing chambers per disc	12		12 + 12	
Power consumption (kW) Without L.F. – With L.F.	12,7 – 16,7			
Compressed air for powder ejection (at 7 bar)	120 NI/min.			

